


Summary

Name	UBP11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P51784
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LCFRFRDQNP EVAVEGRLPI SHSCVGCRRR RTAMATVAAN PAAAAA AVAA AAAVTEDREP QHEELPGLDS QWRQIENGES GRERPLRAGE SWFLVEKHWHY KQWEAYVQGG DQDSSTFPGC INNATLFQDE INWRLKEGLV EGEDYVLLPA AAWHYLVSWY GLEHGQPPIE RKVIELPNIQ KVEVYPVELL LVRHNDLGKS HTVQFSHTDS IGLVLR TARE RFLVEPQEDT RLWAKNSEGS LDRLYDTHIT VLDAALETGQ LIIMETRKKD GTWPSAQLHV MNNNMSEEDE

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 11 (EC 3.4.19.12) (Deubiquitinating enzyme 11) (Ubiquitin thioesterase 11) (Ubiquitin-specific-processing protease 11)
--------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------

Product Name: Recombinant Human USP11 Protein
Catalog #: PEH8631



Note

For Research Use Only , Not for Diagnostic Use