


Product Name: Recombinant Human UBE4A Protein
Catalog #: PEH8633



Summary

Name	UBE4A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14139
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DELPASPDDS DNSVSESLDE FDYSVAEISR SFRSQQEICE QLNINHMIQR IFLITLDNSD PSLKSGNGIP SRCVYLEEMA VELEDQDWLD MSNVEQALFA RLLQDPGNH LINMTSSTTL NLSADRDAGE RHIFCYLYSC FQRAKEEITK VPENLLPFAV QCRNLTVSNT RTVLLTPEIY VDQNIHEQLV DLMLEAIQGA HFEDVTEFLE EVIEALILDE EVRTFPEVMI PVFDILLGRI KDLELCQILL YAYLDILLYF TRQKDMAKVF VEYIQPKDPT NGQMYQKTLL GVILSISCLL KTPGVVENHG YFLNPSRSSP QEIKVQEANI HQFMAQFHEK IYQMLKNLLQ LSPETKHIL SWLGNCLHAN AGRTKIWANQ MPEIFFQMYA SDAFFLNLGA

image



Background

Alternative Names	Ubiquitin conjugation factor E4 A (EC 2.3.2.27) (RING-type E3 ubiquitin
--------------------------	---

Product Name: Recombinant Human UBE4A Protein
Catalog #: PEH8633



transferase E4 A)

Note

For Research Use Only , Not for Diagnostic Use