


Product Name: Recombinant Human VPS41 Protein
Catalog #: PEH8635



Summary

Name	VPS41_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49754
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLDNEDKISI KKVVEELED R PELQHVVYLHK LFKRDHHKGQ RYHEKQISLY AEYDRPNLLP FLRDSTHCPL EKALEICQQR NFVEETVYLL SRMGNSRSAL KMIMEELHDV DKAIEFAKEQ DDGELWEDLI

image



Background

Alternative Names Vacuolar protein sorting-associated protein 41 homolog (S53)

Note

For Research Use Only , Not for Diagnostic Use