


Product Name: Recombinant Human ZPR1 Protein
Catalog #: PEH8637



Summary

Name	ZPR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75312
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VRYTLSVRAL EDMNREVVKT DSAATRIPEL DFEIPAFSQK GALTIVEGLI TRAISSGLEQD QPARRANKDA TAERIDEFIV KLKELKQVAS PFTLIIDDP GNSFVENPHA PQKDDALVIT HYNRTRQQEE MLGLQEEAPA EKPEEEDLRN EVLQFSTNCP ECNAPAQTNM KLVQIPHFKV VIIMATNCEN CGHRTNEVKS GGAVEPLGTR ITLHITDASD MTRDLLKSET CSVEIPELEF ELGMAVLGGK FTTLEGLLKD IRELVTKNPF TLGDSSNPGQ TERLQEFSSQK MDQIIEGNMK AHFIMDDPAG NSYLQNVYAP EDDPEMKVER

image



Background

Alternative Names Zinc finger protein ZPR1 (Zinc finger protein 259)

Product Name: Recombinant Human ZPR1 Protein
Catalog #: PEH8637



Note

For Research Use Only , Not for Diagnostic Use