


Product Name: Recombinant Human YTHDC2 Protein
Catalog #: PEH8641



Summary

Name	YTDC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H6S0
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDILRTTGYT NKEMLKYKKE KQEEKQQT LTEWYSAQEN SFKPESQRQR TVLNVTD EYD LLDDGGDAVF SQLTEKDVNC LEPWLIKEMD ACLSDIWLHK DIDAF AQVFH LILTENVSD YRHSETSATA LMVAAGRGFA SQVEQLISMG ANVH SKASNG WMALDWAKHF GQTEIVDLLE SYSATLEFGN LDESSLVQTN GSDLSAEDRE LLKAYHHSFD

image



Background

Alternative Names	3'-5' RNA helicase YTHDC2 (EC 3.6.4.13) (YTH domain-containing protein 2) (hYTHDC2)
--------------------------	-------------------------------------------------------------------------------------

Note

Product Name: Recombinant Human YTHDC2 Protein
Catalog #: PEH8641



For Research Use Only , Not for Diagnostic Use