


Product Name: Recombinant Human ACTR3B Protein
Catalog #: PEH8644



Summary

Name	ARP3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P1U1
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YAGNTEPQFI IPSCIAIRES AKVVDQAQRR VLRGVDDLDF FIGDEAIDKP TYATKWPIRH GIIEDWDLME RFMEQVVKY LRAEPEDHYF LMTEPPLNTP ENREYLAEIM FESFNVPGLY IAVQAVLALA ASWTSRQVGE RTLTGIVIDS GDGVTHVIPV AEGYVIGSCI KHIPIAGRDI TYFIQQLLRE REVGIPPEQS LETAKAIKEK YCYICPDIVK EFAKYDVDPR KWIKQYTGIN AINQKKFVID VGYERFLGPE IFFHPEFANP DMESISDVV DEVIQNCPID VRRPLYKNV LSGGSTMFRD FGRRQLRDLK RVVDARLRLS EELSGGRIKP KPVEVQVWTH HMQRVAVWFG GSMLASTPEF FQVCHTKKDY EEYGPSICRH

image



Background

Alternative Names Actin-related protein 3B (ARP3-beta) (Actin-like protein 3B) (Actin-related

Product Name: Recombinant Human ACTR3B Protein
Catalog #: PEH8644



protein ARP4)

Note

For Research Use Only , Not for Diagnostic Use