


Product Name: Recombinant Human ACSM2A Protein
Catalog #: PEH8646



Summary

Name	ACS2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q08AH3
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LDHWADMEKA GKRLPSPALW WVNGKGKELM WNFRELSNS QQAANVLSGA CGLQRGDRVA VVLPRVPEWW LVILGCIRAG LIFMPGTIQM KSTDILYRLQ MSKAKAIVAG DEVIQEVDTV ASECPSLRIK LLVSEKSCDG WLNFKKLLNE ASTTHHCVET GSQEASAIYF TSGTSGLPKM AEHSYSSLGL KAKMDAGWTG LQASDIMWTI SDTGWILNIL CSLMEPWALG ACTFVHLLPK FDPLVILKTL SSYPIKSMMG APIVYRMLLQ QDLSSYKFPH LQNCVTVGES LLPETLENWR AQTGLDIRES YGQTETGLTC MVSKTMKIKP GYMGTAASCY DVQIIDDKGN VLPPGTEGDI GIRVKPIRPI GIFSGYVDNP DKTAANIRGD FWLLGDRGIK DEDGYFQFMG RANDIINSSG YRIGPSEVEN ALMEHPAVVE TAVISSPDPV RGEVVKAFVV LASQFLSHDP EQLTKELQQH VKSVTAPYKY PRKIEFVLNL PKTVTGKIQR

image



Background

Product Name: Recombinant Human ACSM2A Protein
Catalog #: PEH8646



Alternative Names

Acyl-coenzyme A synthetase ACSM2A, mitochondrial (EC 6.2.1.2) (Acyl-CoA synthetase medium-chain family member 2A) (Benzoate--CoA ligase) (EC 6.2.1.25) (Butyrate--CoA ligase 2A) (Butyryl-coenzyme A synthetase 2A) (Middle-chain acyl-CoA synthetase 2A)

Note

For Research Use Only , Not for Diagnostic Use