


Product Name: Recombinant Human ACSM2B Protein
Catalog #: PEH8647



Summary

Name	ACS2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q68CK6
Host	E.coli
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LDHWADMEKA GKRLPSPALW WVNGKGKELM WNFRELSSENS QQAANILSGA CGLQRGDRVA VMLPRVPEWW LVILGCIRAG LIFMPGTIQM KSTDILYRLQ MSKAKAIVAG DEVIQEVDTV ASECPSLRIK LLVSEKSCDG WLNFKKLLNE ASTTHHCVET GSQEASAIYF TSGTSGLPKM AEHSYSSLGL KAKMDAGWTG LQASDIMWTI SDTGWILNIL GSLLESWTLG ACTFVHLLPK FDPLVILKTL SSYPIKSMMG APIVYRMLLQ QDLSSYKFPH LQNCLAGGES LLPETLENWR AQTGLDIREF YGQTETGLTC MVSKTMKIKP GYMGTAAASY DVQVIDDKGN VLPPGTEGDI GIRVKPIRPI GIFSGYVENP DKTAANIRGD FWLLGDRGIK DEDGYFQFMG RADDIINSSG YRIGPSEVEN ALMKHPAVVE TAVISSPDPV RGEVVKAFVI LASQFLSHDP EQLTKELQQH VKSVTAPYKY PRKIEFVLNL

image



Background

Product Name: Recombinant Human ACSM2B Protein
Catalog #: PEH8647



Alternative Names

Acyl-coenzyme A synthetase ACSM2B, mitochondrial (EC 6.2.1.2) (Acyl-CoA synthetase medium-chain family member 2B) (Benzoate--CoA ligase) (EC 6.2.1.25) (Butyrate--CoA ligase 2B) (Butyryl-coenzyme A synthetase 2B) (Middle-chain acyl-CoA synthetase 2B) (Xenobiotic/medium-chain fatty acid-CoA ligase HXM-A)

Note

For Research Use Only , Not for Diagnostic Use