


## Summary

<b>Name</b>	ADAL_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q6DHV7
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	36 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	KTDFYSELPK VELHAHLNGS ISSHTMKKLI AQKPDLEIHD QMTVIDKGKK RTLEECFQMF QTIHQLTSSP EDILMVTKDV IKEFADDGVK YLELRSTPRR ENATGMTKKT YVESILEGIK QSKQENLDID VRYLIAVDRR GGPLVAKETV KLAEEFFLST EGTVLGLDLS GDPTVGQAKD FLEPLLEAKK AGLKLALHLS EIPNQKQKQ ILLDLLPDRI GHGTFLNSGE GGSLDLVDFV RQHRIPLELC LTSNVKSQTV PSYDQHFGF WYSIAHPSVI CTDDKGVFAT HLSQEYQLAA

## image



## Background

<b>Alternative Names</b>	Adenosine deaminase-like protein (EC 3.5.4.-) (Adenosine deaminase-like protein isoform 1) (N6-mAMP deaminase) (HsMAPDA) (N6-methyl-AMP aminohydrolase)
--------------------------	---

**Product Name: Recombinant Human ADAL Protein**  
**Catalog #: PEH8649**

---



**Note**

For Research Use Only , Not for Diagnostic Use