


Product Name: Recombinant Human AKR7L Protein
Catalog #: PEH8651



Summary

Name	ARK74_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NHP1
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ATVLGAMEMG RRMDAPTSAA VTRAFLEGRH TEIDTAFLYS DGQSETILGG LGLRMGSSDC RVKIATKANP WIGNSLKPDS VRSQLETSLK RLQCPRVDLF YLHAPDHSAP VEETLRACHQ LHQEGKFVEL GLSNYAWEV AEICTLCKSN GWILPTVYQG MYSATTRQVE TELFPCLRHF GLRFYAYNPL AGGLLTGKYK YEDKDGKQPV GRFFGTQWAE IYRNHFWEKH HFEGIALVEK ALQAAYGASA PSMTSAALRW MYHHSQLQGA HGDAVILGMS SLEQLEQNLA AAEEGPLLEPA VVDAFNQAWH

image



Background

Alternative Names	Aflatoxin B1 aldehyde reductase member 4 (EC 1.-.-.-) (AFB1 aldehyde reductase 3) (AFB1-AR 3) (Aldoketoreductase 7-like)
--------------------------	--

Product Name: Recombinant Human AKR7L Protein
Catalog #: PEH8651



Note

For Research Use Only , Not for Diagnostic Use