


**Product Name: Recombinant Human AFTPH Protein**  
**Catalog #: PEH8652**



## Summary

<b>Name</b>	AFTIN_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q6ULP2
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	45 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	FSENVDSLTS FSIKNGNDK DITAEISAPV KGQSDVLLST TSKEIISSEM LATSIDGMER PGNLNKVVVEQ RQNVGTLESF SPGDFRTNMN VVHQNKQLES CNGEKPPCLE ILTNGFAVLE TVNPQGTDDL DNVADSKGRK PLSTHSTEYN LDSVPSPAEE FADFATFSKK ERIQLEEIEC AVLNDREALT IRENNKINRV NELNSVKEVA LGRSLDNKGD TDGEDQVCVS EISIVTNRGF SVEKQGLPTL QQDEFLLQSGV QSKAWSLVDS ADNSEAIRRE QCKTEEKLDL LTSKCAHLCM DSVKTSDDDEV GSPKEESRKF TNFQSPNIDP TEENDLDDSL SVKNGDSSND FVTCNDINED DFGDFGDFGS ASGSTPPFVT GTQDSMSDAT

## image



## Background

**Alternative Names** Aftiphilin

**Product Name: Recombinant Human AFTPH Protein**  
**Catalog #: PEH8652**

---



**Note**

For Research Use Only , Not for Diagnostic Use