


Product Name: Recombinant Human AKR1B15 Protein
Catalog #: PEH8653



Summary

Name	AK1BF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	C9JRZ8
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKMPIVGLGT WRSLLGKVKE AVKVAIDAEY RHIDCAYFYE NQHEVGAEIQ EKIQEKAVMR EDLFIVSKVW PTFFERPLVR KAFEKTLKDL KLSYLDVYLI HWPQGFKTGD DFFPKDDKGN MISGKGTFDL AWEAMEELVD EGLVKALGVS NFNHFQIERL LNKPGKLYKP VTNQVECHPY LTQEKLQYC HSKGITVTAY SPLGSPDRPW AKPEDPSLLE DPKIKEIAAK HKKTTAQVLI RFHIQRNVTV IPKSMTPAHI VENIQVDFDK LSDEEMATIL SFNRNWRAFD KKEFSHLEDF

image



Background

Alternative Names	Aldo-keto reductase family 1 member B15 (EC 1.1.1.-) (EC 1.1.1.300) (EC 1.1.1.54) (Estradiol 17-beta-dehydrogenase AKR1B15) (Farnesol dehydrogenase) (EC 1.1.1.216) (Testosterone 17beta-dehydrogenase) (EC
--------------------------	---

Product Name: Recombinant Human AKR1B15 Protein
Catalog #: PEH8653



1.1.1.64)

Note

For Research Use Only , Not for Diagnostic Use