

Product Name: Recombinant Human ZFAND2B Protein
Catalog #: PEH8656



Summary

Name	ZFN2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WV99
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSRAGLAAIS RAQAVASTST VPSPSQTMPST CTSPSRATTR SPSWTAPPVI ALQNGLSEDE ALQRALEMSL AETKPQVPSC QEEEDLALAQ ALSASEAEYQ RQQAQSRSSK

image



Background

Alternative Names AN1-type zinc finger protein 2B (Arsenite-inducible RNA-associated protein-like protein) (AIRAP-like protein)

Note

For Research Use Only , Not for Diagnostic Use