


Product Name: Recombinant Human ANKHD1 Protein
Catalog #: PEH8660



Summary

Name	ANKH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWZ3
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPDRTSQENS PALLGVQKGT SKQKSSSLQV ADQDLLPSFH PYQPLECIVE ETEGKLNELG QRISAIEKAQ LKSLELIQGE PLNKDKIEEL KKNREEQVQK KKKILKELQK VERQLQMKTQ QQFTKEYLET KGQKDTVSLH QQCSHRGVFP EGEGDGS LPE DHFSELPQVD TILFKDNDVD DEQQSPPSAE QIDFVPVQPL SSPQCNFSSD

image



Background

Alternative Names Ankyrin repeat and KH domain-containing protein 1 (HIV-1 Vpr-binding ankyrin repeat protein) (Multiple ankyrin repeats single KH domain) (hMASK)

Note

Product Name: Recombinant Human ANKHD1 Protein
Catalog #: PEH8660



For Research Use Only , Not for Diagnostic Use