


Product Name: Recombinant Human ANKS1A Protein
Catalog #: PEH8661



Summary

Name	ANS1A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92625
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LTALDTVREL PSQKSQQIAA LIEDHMTGKR STKEVDKTPP PQPLISSMD SISQKSQGDV EKAVTELIID FDANAEIEGP YEALYNAISC HSLDSMASGR SSDQDSTNKE AEAAGVKPAG VRPRERPPPP AKPPPDEEEE DHIDKKYFPL TASEVLSMRP RIHGSAAREE DEHPYELLT AETKKVVLVD GKTKDHRSS SSRSQDSAEG QDGQVPEQFS GLLHGSSPVC EVGQDPFQLL CTAGQSHPDG SPQQGACHKA

image



Background

Alternative Names Ankyrin repeat and SAM domain-containing protein 1A (Odin)

Product Name: Recombinant Human ANKS1A Protein
Catalog #: PEH8661



Note

For Research Use Only , Not for Diagnostic Use