


Product Name: Recombinant Human ASB3 Protein
Catalog #: PEH8662



Summary

Name	ASB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y575
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CSTVGLAARE GNVKVLRLKLL KKGSRVDVAD NRGWMPHEA AYHNSVECLQ MLINADSEN YIKMKTFFEGF CALHLAASQG HWKIVQILLE AGADPNATTL EETPLFLAV ENGQIDVLRLL LQHGAVNG SHSMCGWNSL HQASFQENAE IIKLLLRKGA NKECQDDFGI TPLFVAAQYG KLESLSILIS SGANVNCQAL DKATPLFIAA QEGHTKCVEL

image



Background

Alternative Names Ankyrin repeat and SOCS box protein 3 (ASB-3)

Note

Product Name: Recombinant Human ASB3 Protein
Catalog #: PEH8662



For Research Use Only , Not for Diagnostic Use