


Product Name: Recombinant Human ANKZF1 Protein
Catalog #: PEH8664



Summary

Name	ANKZ1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	Q9H8Y5
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVLPKRRRRK RNKKEKSRDQ EAGAHRTLLQ QTQEEEPSTQ SSQAVAAPLG PLLDEAKAPG QPELWNALLA ACRAGDVGVL KLQLAPSPAD PRVLSLLSAP LGSGGFTLLH AAAAAGRGSV VRLLEAGAD PTVQDSRARP PYTVAADKST RNEFRRFMEK NPDAYDYNKA QVPGPLTPEM EARQATRKRE QKAARRQREE QQQRQQEQEE REREEQRRFA ALSDREKRAL AAERRLAAQL GAPTSPIPDS AIVNTRRCWS CGASLQGLTP FHYLDFSFCs TRCLQDHRRO

image



Background

Alternative Names Ankyrin repeat and zinc finger domain-containing protein 1 (Zinc finger protein 744)

Product Name: Recombinant Human ANKZF1 Protein
Catalog #: PEH8664



Note

For Research Use Only , Not for Diagnostic Use