


**Product Name: Recombinant Human ANKRD50 Protein**  
**Catalog #: PEH8665**



## Summary

<b>Name</b>	ANR50_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-His
<b>Accession #</b>	Q9ULJ7
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	37 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GASSLNGCSP SPVHTMEQKP LQSLSSKVQS LTIKSNSSGS TGGGDMQPSL RGLPNGPTHA FSSPESPDS TVDRQKSSLS NNSLKSSKNS SLRTTSSTAT AQTVPIDSFH NLSFTEQIQQ HSLPRRSRQ SIVSPSSTTQ SLGQSHNSPS SEFEWSQVKP SLKSTKASKG GKSENSAKSG SAGKKAKQSN SSQPKVLEYE MTQFDRRGPI AKSGTAAPPK QMPAESQCKI MIPSAAQQEIG RSQQQFLIHQ QSGEQKKRNG IMTNPNYHLQ SNQVFLGRVS VPRTMQDRGH QEVLEGYPSS ETELSLKQAL KLQIEGSDPS

## image



SDS

## Background

**Alternative Names** Ankyrin repeat domain-containing protein 50

**Product Name: Recombinant Human ANKRD50 Protein**  
**Catalog #: PEH8665**



---

**Note**

For Research Use Only , Not for Diagnostic Use