


Product Name: Recombinant Human ANKRD6 Protein
Catalog #: PEH8666



Summary

Name	ANKR6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9Y2G4
Host	E.coli
Species	Human
Molecular Mass	72 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QTPLETARYH NNPEVALLT KAPQVLRFSR GRSLRKKRER LKEERRAQSV PRDEVAQSKG SVSAGDTPSS EQAVARKEEA REEFLSASPE PRAKDDRRRK SRPKVSAFSD PTPPADQQPG HQKNLHAHNH PKKRNRRHCS SPPPPHEFRA YQLYTLYRGK DGKVMQAPIN GCRCEPLINK LENQLEATVE EIKAELGSVQ DKMNTKLGQM ENKTQHQMRV LDKLMVERLS AERTECLNRL QQHSDTEKHE

image



Background

Alternative Names Ankyrin repeat domain-containing protein 6 (Diversin)

Note

Product Name: Recombinant Human ANKRD6 Protein
Catalog #: PEH8666



For Research Use Only , Not for Diagnostic Use