

Product Name: Recombinant Human AP1S1 Protein
Catalog #: PEH8668



Summary

Name	AP1S1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P61966
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IEGQDNELIT LELIHRYVEL LDKYFGSVCE LDIIFNFEKA YFILDEFLMG GDVQDTSKKS VLKAIQADL LQEEDESPRS

image



Background

Alternative Names	AP-1 complex subunit sigma-1A (Adaptor protein complex AP-1 subunit sigma-1A) (Adaptor-related protein complex 1 subunit sigma-1A) (Clathrin assembly protein complex 1 sigma-1A small chain) (Clathrin coat assembly protein AP19) (Golgi adaptor HA1/AP1 adaptin sigma-1A subunit) (HA1 19 kDa subunit) (Sigma 1a subunit of AP-1 clathrin) (Sigma-adaptin 1A) (Sigma1A-adaptin)
--------------------------	--

Product Name: Recombinant Human AP1S1 Protein
Catalog #: PEH8668



Note

For Research Use Only , Not for Diagnostic Use