


Product Name: Recombinant Human AP3B2 Protein
Catalog #: PEH8669



Summary

Name	AP3B2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13367
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DQLIEVIEKL LADKTTLVAG SVVMAFEEVC PERIDLIHKN YRKLCNLLID VEEWGQVVII SMLTRYARTQ FLSPTQNESL LEENAEKAFY GSEEDEAKGA GSEETAAAAA

image



Background

Alternative Names	AP-3 complex subunit beta-2 (Adaptor protein complex AP-3 subunit beta-2) (Adaptor-related protein complex 3 subunit beta-2) (Beta-3B-adaptin) (Clathrin assembly protein complex 3 beta-2 large chain) (Neuron-specific vesicle coat protein beta-NAP)
--------------------------	---

Note

Product Name: Recombinant Human AP3B2 Protein
Catalog #: PEH8669



For Research Use Only , Not for Diagnostic Use