


Product Name: Recombinant Human APOL5 Protein
Catalog #: PEH8670



Summary

Name	APOL5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BWW9
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPGSKVLPGL GEGCKEMWLR KVIYGGGEVWG KSPEPEFPSL VNLCQSWKIN NLMSTVHSDE AGMLSYFLFE ELMRCDKDSM PDGNLSEEEK LFLSYFPLHK FELEQNIKEL NTLADQVDTT HELLTKTSLV ASSSGAVSGV MNILGLALAP VTAGGSLMLS ATGTGLGAAA AITNIVTNVL ENRSNSAARD KASRLGPLTT SHEAFGGINW SEIEAAGFCV NKCVKAIQGI KDLHAYQMAK SNSGFMAMVK NFVAKRHIPF WTARGVQRAF EGTTLAMTNG AWVMGAAGAG FLLMKDMSSF LQSWKHLEDG ARTETAELR ALAKKLEQEL DRLTQHHRHL POKASQTCSS

image



Background

Alternative Names Apolipoprotein L5 (Apolipoprotein L-V) (ApoL-V)

Product Name: Recombinant Human APOL5 Protein
Catalog #: PEH8670



Note

For Research Use Only , Not for Diagnostic Use