


Product Name: Recombinant Human ASAP3 Protein
Catalog #: PEH8672



Summary

Name	ASAP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TDY4
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSVLEAFANG QDFGQPLPGP DAQAPEELVL HLAVKVANQA SLPLVDFIIQ NGGHLDAKAA DGNTALHYAA LYNQPDCLKL LLKGRALVGT VNEAGETALD IARKKHHKEC EELLEQAQAG TFAFPLHVDY SWVISTEPGS DSEEDDEEEKR

image



Background

Alternative Names	Arf-GAP with SH3 domain, ANK repeat and PH domain-containing protein 3 (Development and differentiation-enhancing factor-like 1) (Protein up-regulated in liver cancer 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ASAP3 Protein
Catalog #: PEH8672

