


Product Name: Recombinant Human DDX24 Protein
Catalog #: PEH8676



Summary

Name	DDX24_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9GZR7
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGIKVVGKWK EVKIDPNMFA DGQMDDLVCF EELTDYQLVS PAKNPSSLFS KEAPKRKAQA VSEEEEEEG KSSSPKKKIK LKKSKNVATE GTSTQKEFEV KDPELEAQGD DMVCDDPEAG EMTSENLVQT APKKKKNKGG

image



Background

Alternative Names ATP-dependent RNA helicase DDX24 (EC 3.6.4.13) (DEAD box protein 24)

Note

For Research Use Only , Not for Diagnostic Use