


Product Name: Recombinant Human BAG2 Protein
Catalog #: PEH8679



Summary

Name	BAG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95816
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NEGRFCRSSS MADRSSRLE SLDQLELRVE ALREAATAVE QEKEILLEMI HSIQNSQDMR QISDGEREEL NLTANRLMGR TLTVEVSVET IRNPQQQESL KHATRIIDEV VNKFLDDLGN AKSHLMSLYS ACSSEVPHGP VDQKFQSIVI GCALEDAQKI KRRLETLRN IENSDKAIKL LEHSGAGASK

image



Background

Alternative Names BAG family molecular chaperone regulator 2 (BAG-2) (Bcl-2-associated athanogene 2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BAG2 Protein
Catalog #: PEH8679

