



Product Name: Recombinant Human EPB41L4A Protein
Catalog #: PEH8680



Summary

Name	E41LA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HCS5
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VKSFKKAKNE NSPDTQRSKS HAPWEENGPO SGLYNPSDR TKSPKFPYTR RRNPSCGSDN DSVQPVRRRK AHNSGEDSDL KQRRRSR SRC NTSSGSESEN SNREYRKKRN RIRQENDMVD SAPQWEAVLR RQKEKNQADP NNRRSRHRSR SRSPDIQAKE ELWKHIQKEL VDPSGLSEEQ LKEIPYTKIE TQGDPIRIRH SHSPRSYRQY RRSQCSDGER SVLSEVNSKT DLVPPLPVTR SSDAQGSGDA TVHQRRNGSK DSLMECKPQT STNNLAGKHT AKTIKTIQAS

image



Background

Alternative Names	Band 4.1-like protein 4A (Erythrocyte membrane protein band 4.1-like 4A) (Protein NBL4)
--------------------------	---

Product Name: Recombinant Human EPB41L4A Protein
Catalog #: PEH8680



Note

For Research Use Only , Not for Diagnostic Use