


Product Name: Recombinant Human BEGAIN Protein
Catalog #: PEH8690



Summary

Name	BEGIN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BUH8
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGELRKRLSY TTHKLEKLET EFDSTRHYLE IELRRAQEEL EKVTEKLRLRI QSNYMALQRI NQELEDKLYR MGQHYEEER ALSHEIVALN SHLLEAKVTI DKLSEDNELY RKDCNLAAQL LQCSQTYGRV HKVSELP SDF QERVSLHMEK HGCSLP S PLC HPAYADSVPT CVIAKVLEKP DPASLSSRLS DASARDLAFC DGVEKPGPRP PYKGDYICSD TALYCPEERR RDRRPSVDAP VTDVGFLRAQ

image



Background

Alternative Names Brain-enriched guanylate kinase-associated protein

Note

Product Name: Recombinant Human BEGAIN Protein
Catalog #: PEH8690



For Research Use Only , Not for Diagnostic Use