


Summary

| | |
|------------------------------|---|
| Name | B12L2_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q6UXY1 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 29 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | GSGSYGTEPD ARPASQLEPD RRSLPRTPSA SSLYSGSAQS SRSNSFGERP GGGGGARRVR ALVSHSEGAN HTLLRFSGAGD VVEVLVPEAQ NGWLYGKLEG SSASGWFPEA YVKALEEGPV NPMTPTVTPMT SMTSMSPMTP MNPGNELPSR SYPLRGSHSL DDLLDRPGNS IAPSEYWDGQ SRSRTPSRVP SRAPSPAPP LPSSRRSSMG STAVATDVKK LMSSEQYPPQ ELFPRGTNPF ATVKLRPTIT |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 2 (BAI1-associated protein 2-like protein 2) (Planar intestinal- and kidney-specific BAR domain protein) (Pinkbar) |
|--------------------------|--|

Product Name: Recombinant Human BAIAP2L2 Protein
Catalog #: PEH8691



Note

For Research Use Only , Not for Diagnostic Use