


Product Name: Recombinant Human BICRAL Protein
Catalog #: PEH8693



Summary

Name	BICRL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6AI39
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GNISKKTECL GRALKFDKVG LVQYQSTSEE KASRREPLKA SQCSPGPEGH RKTSSRSDHG TESKLSSILA DSHLEMTCNN SFQDKSLRNS PKNEVLHTDI MKGSGEPQPD LQLTKSLETT FKNILELKKA GRQPQSDPTV SGSVELDFPN FSPMASQENC LEKFIPDHSE

image



Background

Alternative Names	BRD4-interacting chromatin-remodeling complex-associated protein-like (Glioma tumor suppressor candidate region gene 1 protein-like)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BICRAL Protein
Catalog #: PEH8693

