

Product Name: Recombinant Human KCTD2 Protein
Catalog #: PEH8697



Summary

Name	KCTD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14681
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGLGGGGGSG VGDGGGPVRG PPSRPAGPT PRGHGRPAAA VAQPLEPGPG PPERAGGGGA ARWVRLNVGG TYFVTTRQTL GREPKSFLCR LCCQEDPELD SDKDETGAYL IDRPTYFGP ILNYLRHGKL IITKELAEEG VLEEAEFYNI ASLVRLVKER IRDNENRTSQ GPVKHVYRVL QCQEEELTQM VSTMSDGWKF EQLISIGSSY NYGNEDQAEF LCVVSRELNN STNGIVIEPS EKAKILQERG

image



Background

Alternative Names	BTB/POZ domain-containing protein KCTD2 (Potassium channel tetramerization domain-containing protein 2)
--------------------------	---

Note

Product Name: Recombinant Human KCTD2 Protein
Catalog #: PEH8697



For Research Use Only , Not for Diagnostic Use