


Product Name: Recombinant Human CABP2 Protein
Catalog #: PEH8699



Summary

Name	CABP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NPB3
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGNCAKRPWR RGPKDPLQWL GSPPRGSCPS PSSSPKEQGD PAPGVQGYSV LNSLVGPACI FLRPSIAATQ LDRELRP EELQVAFQEF DRDRDGYIGC RELGACMRTL GYMPTMELI EISQQISGGK VDFEDFVELM GPKLLAETAD MIGVREL RDA FREFDTNGDG RISVGELRAA LKALLGERLS QREVDEILQD VDLNGDGLVD FEEFVRMMSR

image



Background

Alternative Names Calcium-binding protein 2 (CaBP2)

Note

Product Name: Recombinant Human CABP2 Protein
Catalog #: PEH8699



For Research Use Only , Not for Diagnostic Use