


Product Name: Recombinant Human CLIP3 Protein
Catalog #: PEH8705



Summary

Name	CLIP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96DZ5
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLDGQKTGTL RFCGTTEFAS GQWVGVELDE PEGKNDGSVG GVRVFICPPK QGLFASVSKI SKAVDAPPSS VTSTPRTPRM DFSRVTKGGR REHKGKKKTP SSPSLGSLQQ RDGAKAEVGD

image



Background

Alternative Names	CAP-Gly domain-containing linker protein 3 (Cytoplasmic linker protein 170-related 59 kDa protein) (CLIP-170-related 59 kDa protein) (CLIPR-59)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use