


Product Name: Recombinant Human CASKIN1 Protein
Catalog #: PEH8706



Summary

Name	CSK11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WXD9
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPKKLEPASE GKSSEAVSQW LTAFLQLQLYA PNFISAGYDL PTISRMTPEL LTAIGVTKPG HRKKIAAEIS GLSIPDWLPE HKPANLAVWL SMIGLAQYYK VLVDNGYENI DFITDITWED LQEIGITKLG HQKKLMLAVR KLAELQKAAY AKYEGGPLRR KAPQSLEVMA IESPPPPEPT PADCCSPKMT TFQDSELSDE LQAAMTGPAE

image



Background

Alternative Names Caskin-1 (CASK-interacting protein 1)

Note

Product Name: Recombinant Human CASKIN1 Protein
Catalog #: PEH8706



For Research Use Only , Not for Diagnostic Use