


Summary

| | |
|------------------------------|--|
| Name | CECR2_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-His |
| Accession # | Q9BXF3 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 32 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | RAKRRKLREE RAWLLAQGKE LPPELSHLDP NSPMREEKKT KDLFELDDDF TAMYKVLDDV KAHKDSWPFL EPVDESYAPN YYQIIKAPMD ISSMEKKLNG GLYCTKEEFV NDMKTMFRNC RKYNGESSEY TKMSDNLERC FHRAMMKHFP GEDGDTDEEF WIREDEKREK RRSRAGRSGG SHVWTRSRDP EGSSRKQQPM ENGGKSLPPT RRAPSSGDDQ SSSSTQPPRE VGTSNRGRGFS HPLHCGGTPS QAPFLNQMRP AVPGTFGPLR |

image



Background

Alternative Names Cat eye syndrome critical region protein 2

Product Name: Recombinant Human CECR2 Protein
Catalog #: PEH8708



Note

For Research Use Only , Not for Diagnostic Use