


Product Name: Recombinant Human CATSPERZ Protein
Catalog #: PEH8709



Summary

Name	CTSRZ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9NTU4
Host	E.coli
Species	Human
Molecular Mass	65 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEEKPSKVSL KSSDRQGSDE ESVHSDTRDL WTTTTLSQAQ LNMPLSEVCE GFDEEGRNIS KTRGWHPGR GSLDEGYKAS HKPEELDEHA LVELELHRGS SMEINLGEKD TASQIEAEKS SSMSLNIK HMPHRAYWAE QQSRLPLPLM ELMENEALEI LTKALRSYQL GIGRDHFLTK ELQRYIEGLK KRRSKRLVYN

image



Background

Alternative Names Cation channel sperm-associated protein subunit zeta (CatSper-zeta) (CatSperzeta) (Testis-expressed protein 40)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CATSPERZ Protein
Catalog #: PEH8709

