



Summary

Name	CNTLN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NXG0
Host	E.coli
Species	Human
Molecular Mass	62 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GAPLPLPQES DPDYSAQVPH RPSLSSLETL MVSQKSEIEY LQEKLKIANE KLSNISANK GFSRKSIMTS AEGKHKEPPV KRSRSLSPKS SFTDSEELQK LRKAERKIEN LEKALQLKSQ ENDELDAHE KRKERLQMLQ TNYRAVKEQL KQWEEGSGMT EIRKIKRADP QQLRQEDSDA VWNELAYFKR ENQELMIQKM NLEEEDELK VHISIDKAAI QELNRCVAER REEQLFRSGE DDEVKRSTPE KNGKEMLEQT LQKVTELENR LKSF EKRSRK LKEGNKMLK ENDFLKSLLK QQQEDTETRE KELEQIIGS KDVEKENTEL QVKISELETE VTSLRRQVAE ANALRNENEE LINPMEKSHQ SADRAKSEMA TMKVRSGRYD CKTTMTKVKF KAAKNC SVG RHHTVLNHSI KVMSNVFENL SKDGWEDVSE SSSDSEAQTS QTLGTIIVET SQKISPTEDG KDQKESDPT E DSQTQGKEIV QTYLNIIDGKT PKDYFHDKNA KKPTFQKKNC

image



Background

Product Name: Recombinant Human CNTLN Protein
Catalog #: PEH8710



Alternative Names Centlein (Centrosomal protein)

Note

For Research Use Only , Not for Diagnostic Use