


Product Name: Recombinant Human CCSAP Protein
Catalog #: PEH8711



Summary

Name	CCSAP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6IQ19
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SASSESSGAG GPAPRCAPPS PPPPVEPATQ EEAERRARGA PEEQDAEAGD AEAEDAEDAA LPALPVKDVE DKPEQQTRTR ETDKSPTSTE PRQQPSALFA RGNRKAVKSP QRSSSKIEN KHPFALYGWG EKQTDGTSQK THNVCASAPV HEIHESALRA KNRRQVEKRK LVAQRQRAHS VDVEKNRKMK ASSENPWMT EYMRCYSARA

image



Background

Alternative Names Centriole, cilia and spindle-associated protein

Note

Product Name: Recombinant Human CCSAP Protein
Catalog #: PEH8711



For Research Use Only , Not for Diagnostic Use