


Product Name: Recombinant Human CNTROB Protein
Catalog #: PEH8712



Summary

Name	CNTRB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N137
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PLGAEDLLSD SSEPPGLNQV SSEVTSQLYA SLRLSRQAEA TARAQLYLPS TSPPEGLDG FAQELSRSL VGLEKNLKKK DGSKHIFEME SVRGQLQTML QTSRDTAYRD PLIPGAGSER REEDSFDS DS TATLLNTRPL QDLSPSSAQ ALEELFPRT SLRPGPLNP PDFQGLRDAL DSEHTRRKHC ERHIQSLQTR VLELQQQLAV AVAADRKKDT MIEQLDKTLA RVVEGWNRHE AERTEVLRGL QEEHQAAELT RSKQQETVTR LEQSLSEAME ALNREQESAR LQQRERETLE EERQALTRL EAEQQRCCVL QEERDAARAG QLSEHRELET LRAALEEERQ TWAQQEHQLK EHYQALQEEQ QAQLEREKEK SQREAQAWE

image



Background

Alternative Names Centrobin (Centrosomal BRCA2-interacting protein) (LYST-interacting

Product Name: Recombinant Human CNTROB Protein
Catalog #: PEH8712



protein 8)

Note

For Research Use Only , Not for Diagnostic Use