


Product Name: Recombinant Human COG5 Protein
Catalog #: PEH8719



Summary

Name	COG5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UP83
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DYLSQGIDLS GIEVIENDLL FIARARLEVE NQAKRLLEQG LETQNPTQVG TALQVFYNLG TLKDTITSVV DGYCATLEEN

image



Background

Alternative Names	Conserved oligomeric Golgi complex subunit 5 (COG complex subunit 5) (13S Golgi transport complex 90 kDa subunit) (GTC-90) (Component of oligomeric Golgi complex 5) (Golgi transport complex 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use