


Product Name: Recombinant Human CORO2A Protein
Catalog #: PEH8720



Summary

Name	COR2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92828
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DRKIRVIDPR AGTVLQEASY KGHRASKVLF LGNLKKLMST GTSRWNNRQV ALWDQDNLSV PLMEEDLDGS SGVLFPPFYDA DTSMLYVVGK GDGNIRYYEV SADKPHLSYL TEYRSYNPQK GIGVMPKRGL DVSSCEIFRF YKLITTKSLI EPISMIVPRR SESYQEDIYP PTAGAQPST AQEWLSGMNR DPILVSLRPG SELLRPHPLP AERPIFNSMA PASPRLNQT

image



Background

Alternative Names Coronin-2A (IR10) (WD repeat-containing protein 2)

Note

Product Name: Recombinant Human CORO2A Protein
Catalog #: PEH8720



For Research Use Only , Not for Diagnostic Use