


Product Name: Recombinant Human COA4 Protein
Catalog #: PEH8723



Summary

Name	COA4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NYJ1
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WTQRVKKDDE EEDPLDQLIS RSGCAASHFA VQECMAQHQD WRQCQPQVQA FKDCMSEQQA RRQEELQRRQ

image



Background

Alternative Names	Cytochrome c oxidase assembly factor 4 homolog, mitochondrial (Coiled-coil-helix-coiled-coil-helix domain-containing protein 8) (E2-induced gene 2 protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use