


Product Name: Recombinant Human PITPNC1 Protein
Catalog #: PEH8727



Summary

Name	PITC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKF7
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TQEATNKKIG IFPPAISISS IPLLSSVRS APSSAPSTPL STDAPEFLSV PKDRPRKKSA PETLTLPDPE KKATLNLPGM HSSDKPCRPK

image



Background

Alternative Names	Cytoplasmic phosphatidylinositol transfer protein 1 (Mammalian rdgB homolog beta) (M-rdgB beta) (MrdgBbeta) (Retinal degeneration B homolog beta) (RdgBbeta)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use