


**Product Name: Recombinant Human DCUN1D2 Protein**  
**Catalog #: PEH8730**



## Summary

<b>Name</b>	DCNL2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q6PH85
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	26 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LHRESMRNAV DKKKLERLYG RYKDPQDENK IGVDGIQQFC DDLSLDPASI SVLVIAWKFR AATQCEFSRK EFLDGMTELG CDSMEKALKAL LPRLEQELKD TAKFKDFYQF TTFKAKNPGQ KGLDLEMAVA YWKLVLSGRF KFLDLWNTFL MEHHKRSIPR DTWNLLDFG NMIADDMSNY DEEGAWPVLII DDFVEYARPV

## image



## Background

**Alternative Names** DCN1-like protein 2 (DCNL2) (DCUN1 domain-containing protein 2) (Defective in cullin neddylation protein 1-like protein 2)

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human DCUN1D2 Protein**  
**Catalog #: PEH8730**

---

