


Product Name: Recombinant Human DCAF5 Protein
Catalog #: PEH8733



Summary

Name	DCAF5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96JK2
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLDGRIEDDS RCLYTHEEYI SLVLNSGSGL SHDYANQSVQ EDPRMMAFFD SLVRREIEGW SSDSDSDLSE STILQLHAGV SERSGYTDSE SSASLPRSP PTVDESADNA FHLGPLRVTT TNTVASTPPT PTCEAASRQ QRLSALRRYQ DKRLLALSNE SDSEENVCEV ELDTDLFPRP

image



Background

Alternative Names	DDB1- and CUL4-associated factor 5 (Breakpoint cluster region protein 2) (BCRP2) (WD repeat-containing protein 22)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DCAF5 Protein
Catalog #: PEH8733

