


Product Name: Recombinant Human GRID2IP Protein
Catalog #: PEH8736



Summary

Name	GRD2I_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A4D2P6
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PQPGPGPICP DSPPSPDPTR PPSRRKLFTF SHPVRSDTD RFLDVLSEQL GPRVTIVDDF LTPENDYEEM SFHDDQGSFV TNERSSASDC ISSSEEGSSL TYSSISDHIP PPPLSPPPPP PLPFHDAKPS

image



Background

Alternative Names Delphilin (Glutamate receptor, ionotropic, delta 2-interacting protein 1)

Note

For Research Use Only , Not for Diagnostic Use