


Product Name: Recombinant Human DPYSL4 Protein
Catalog #: PEH8738



Summary

Name	DPYL4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14531
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ITRWHESIKE ELEALVKEKG VNSFLVMAY KDRQCSDSQ MYEIFSIRD LGALAQVHAE NGDIVEEEQK RLELGITGP EGHVLSHPPEE VEAEAVYRAV TIAKQANCPL YVTKVMSKGA ADAIAQAKRR GVVVFGEPIIT ASLGTDGSHY WSKNWAKAAA FVTSPVNPDP PTTADHLTCL LSSGDLQVTG SAHCTFTTAQ KAVGKDNFAL IPEGTNGIEE RMSMVWEKCV ASGKMDENEF VAVTSTNAAK IFNFYPRKGR VAVGSDADLV IWNPKATKII SAKTHNLNVE YNIFEGVECR GAPAVVISQG RVALEDGKMF

image



Background

Alternative Names Dihydropyrimidinase-related protein 4 (DRP-4) (Collapsin response mediator protein 3) (CRMP-3) (UNC33-like phosphoprotein 4) (ULIP-4)

Product Name: Recombinant Human DPYSL4 Protein
Catalog #: PEH8738



Note

For Research Use Only , Not for Diagnostic Use