


Product Name: Recombinant Human DPH7 Protein
Catalog #: PEH8741



Summary

Name	DPH7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BTV6
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEDRPAGPQN KGGMEVKEPQ VRLGRLFLYS FNDNNSIHPL VEVQRKDTSA ILDMKWCHIP VAGHALLGLA DASGSIQLLR LVESEKSHVL EPLSSLALEE QCLALSLDWS TGKTGRAGDQ PLKIISDST GQLHLLMVNE TRPRLQKVAS WQAHQFEAWI AAFNYWHPEI VYSGGDDGLL RGWDTRVPGK

image



Background

Alternative Names	Diphthine methyltransferase (EC 3.1.1.97) (Diphthamide biosynthesis protein 7) (WD repeat-containing protein 85)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DPH7 Protein
Catalog #: PEH8741

