


Product Name: Recombinant Human DNAH10 Protein
Catalog #: PEH8748



Summary

Name	DYH10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVF4
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KFVEAIEALP LANTPEVFGL HPNAEIGYYT QAARDMWAHL LELQPQTGES SSGISRDDYI GQVAKEIENK MPKVFDLQV RKRLGTGLSP TSVVLLQELE RFNKLVRMT KSLAELQRAL AGEVGMSNEL

image



Background

Alternative Names Dynein heavy chain 10, axonemal (Axonemal beta dynein heavy chain 10)
(Ciliary dynein heavy chain 10)

Note

For Research Use Only , Not for Diagnostic Use