



Product Name: Recombinant Human DRC7 Protein
Catalog #: PEH8753



Summary

Name	DRC7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IY82
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RAQEKRLRE EEERLMEAEK AKPDALHGLR VHSWVLVLSG KREVPENFFI DPFTGHSYST QDEHFLGIES LWNHKNYWIN MQDCWNCKD LIFDLGDPVR WEYMLLGTK SQLSLTEEDD SGINDEDDVE NLGKEDEKKS FDMPSWVEQ IEISPEAFET RCPNGKKVIQ YKRKLEKWA PYLNSNGLVS

image



Background

Alternative Names Dynein regulatory complex subunit 7 (Coiled-coil domain-containing protein 135) (Coiled-coil domain-containing protein lobo homolog)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DRC7 Protein
Catalog #: PEH8753

