


**Product Name: Recombinant Human DTX4 Protein**  
**Catalog #: PEH8754**



## Summary

<b>Name</b>	DTX4_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9Y2E6
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	29 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GVSNTRKTT KKQAKKGKTP EEVLKYLQK VRHPPDEDCT ICMERLTAPS GYKGPQPTVK PDLVGKLSRC GHVYHIYCLV AMYNNGNKDG SLQCPTCKTI YGVKTGTQPP GKMEYHLIPH SLPGHPDCKT IRIYSIPPG IQGPEHPNPG KSFSARGFPR HCYLPDSEKG RKVLKLLLVA WDRRLIFAIG TSSTTGESDT VIWNEVHHKT EFGSNLTGHG YPDANYLDNV LAELAAQGIS

## image



## Background

**Alternative Names** E3 ubiquitin-protein ligase DTX4 (EC 2.3.2.27) (Protein deltex-4) (Deltex4) (RING finger protein 155) (RING-type E3 ubiquitin transferase DTX4)

## Note

**Product Name: Recombinant Human DTX4 Protein**  
**Catalog #: PEH8754**



---

For Research Use Only , Not for Diagnostic Use